

NPi[®]-300 Pupillometer

Modern User Interface

- Simple and intuitive icon-based navigation
- Simplified trending screen displaying NPi and Size trends over 12-hour time windows

Upgraded Keypad

- Easy-to-press buttons designed for fast and simple navigation with or without gloves

NPi-300 Wireless Charging Station

- Completely enclosed wireless charging system with no metal pins or blades for improved reliability and ease of charging
- Durably withstands cleaning with hospital cleaning agents

Incorporated Barcode Scanner

- Instantaneous scanning of 1D or 2D patient barcodes

SmartGuard[®]

- Single patient-use device with RFID memory chip

New On/Off Button

Texturized Plastic Handle

- Ergonomic grips for easy handling



An accurate, reliable and objective system that enhances pupillary assessment to **assist in detecting cerebral insult, guiding treatment and informing prognosis**



prospect
Diagnostics
leaders at the point of care

Prospect Diagnostics Limited

Tel: 01246 292 955 Fax: 01246 292 966
Viking Court, 31 Princess Road, Dronfield S18 2LX
customercare@prospectdiagnostics.co.uk
prospectdiagnostics.co.uk