



Prospect Haemoglobin Analyser (Diaspect™)

Immediate haemoglobin
results
from a hand-held analyser

Results in 1-2 seconds

Reagent-free cuvettes

40 days or 10,000 tests continuous use

2½ years cuvette shelf life from date of production

Fast, easy to use and reliable in all climate conditions

Fast and accurate	<ul style="list-style-type: none"> • Results in 1-2 seconds • Measuring range: 0 - 255 g/L • Imprecision (within run): CV < 1% • Factory calibrated against HiCN reference method according to ICSH • Photometric method with broad spectrum and compensation for turbidity and scattering • No additional calibration required
Easy to use	<ul style="list-style-type: none"> • 'Always on' technology. The analyser is always ready for use • Automatic self-check between every measurement • Long cuvette shelf life: 2.5 years from date of production • Cuvettes unaffected by humidity or temperature and can be used until expiration date after opening the bag
Practical and portable	<ul style="list-style-type: none"> • Sample volume: <10 µL • Capillary or venous blood (EDTA and Heparin anticoagulant can be used) • Maintenance free • Operating ambient temperature: 10 °C to 35 °C • Storage temperature: 0 °C to 50 °C (- 30 °C to 70 °C up to 24 hours maximum.) • Integrated rechargeable lithium-ion battery (40 days / 10,000 tests of continuous use) • Compact: 150 (d) x 90 (w) x 40mm (h) • Lightweight: 180 g • Connect to PC via USB 2.0 or Bluetooth (optional)



1 Collect blood sample



2 Insert the filled cuvette and press down gently



3 Results appear in 1-2 seconds

Distributed by

prospect
Diagnostics Limited
improving clinical outcomes

Prospect Diagnostics Limited
Viking Court, 31 Princess Road
Dronfield, Derbyshire
S18 2LX

☎ 01246 292 955
customer@prospectdiagnostics.co.uk

Manufacturer
DiaSpect Medical GmbH
Von-Cancrin-Str. 1
63877 Sailauf
Germany
☎ +49 6093 996689-60
www.ekfdiagnostics.com

EKF | Diagnostics
for life